

## Patient protection eyewear SOFTCAP

### Product information:

|                         |  |
|-------------------------|--|
| <b>Item no.:</b>        | 600-SOFTCAP-L  |
| <b>Application:</b>     | Diode (561 nm) Diode (630 – 660 nm) Alexandrite (755 nm) Diode (808 nm) ,<br>VLT (Visible light transmission): 0.2 % |
| <b>Filter color:</b>    | green  |
| <b>Filter material:</b> | polycarbonate  |



### Certified protection specifications for EN 207

| Wavelengths | Protection levels |
|-------------|-------------------|
| • 560-820<  | DI LB3 (OD3+)     |
| • 820-840<  | DI LB2 (OD2+) DI  |
| • >840-880< | LB1 (OD1+)        |
|             | PF CE             |

## Properties:

Our SOFTCAPS L (large) patient goggles are made of plastic and are particularly suitable for laser applications with e.g. alexandrite lasers.

Tested according to **EN 207**

- 560 - 820 DI LB3 (OD3+)
- > 820 - 840 DI LB2 (OD2+)
- > 840 - 880 DI LB1 (OD1+)

PF CE