

## Laser/LED disposable eyeshields

### Product information:

**Item no.:** 600-LASP-25-N-TK

**Product weight:** 0 g



### Certified protection specifications for EN 207

#### Wavelengths

- 180-315
- >315-1400
- >1400-11000

#### Protection levels

- D LB7 + IR LB3 + M LB3Y (OD7+)
- D LB6 + IR LB7 + M LB7Y (OD7+)
- DI LB4 (OD4+)

## Properties:

### Patient protection for Laser and LED treatments.

Our latex-free and hypoallergenic eye protection pads provide the patient with light-tight eye protection during treatment. The flat shape simplifies the handling in the treated face areas. The patient's eyebrows and temples remain freely accessible for treatment. The slightly adhesive contact area allows painless positioning and removal. In addition, the flap makes it easier to remove after treatment. The disposable eye protection pads are certified according to European standards.

- 25 pair per pack
- color: white
- also suitable for IPL
- Applied, harmonized norms: EN 207:2017