



IPL protection eyewear ONTOR Shade 2

Product information:

| | |
|-------------------------|------------------------------------|
| Item no.: | 600-K132-ONTO-2-54 |
| Filter color: | green |
| Frame color: | black/white |
| Filter material: | special plastic, scratch resistant |
| Product weight: | 44 g |



Alle Angaben ohne Gewähr.
All information is provided without guarantee.

IPL safety eyewear ONTOR Shade 2 (suitable also for spectacle wearers) offer safety for operator and patient during IPL treatments.

Properties:

Our **IPL safety goggles ONTOR** are lightweight, versatile **safety goggles** made of a plastic frame, which is perfectly suitable as over-glasses for spectacle wearers. The frame offers a very large field of vision and wide side coverage, with excellent wearing comfort due to flexible temples adjustable in length and inclination, and a special soft-rubber padding in the forehead area. **IPL protective goggles ONTOR** provide safety for users and patients during **IPL treatments**.

- Attenuation according to **EN 169: Shade 2**
- suitable for **IPL applications**
- modern sport design
- high wearing acceptance due to perfect fit
- unrestricted field of vision
- suitable for spectacle wearers
- Individually adjustable temple length

- Adjustable inclination
- lightweight (44g)

Included in delivery: microfibre bags with PROTECT-logo, polybag, user's manual, cleaning cloth with PROTECT-logo, glasses head strap