



IPL protection eyewear ONTOR Shade 5 (suitable also for spectacle wearers)

Product information:

Item no.:	600-K135-ONTO-5-54
Filter color:	dark green
Frame color:	black/white
Filter material:	special plastic
Product weight:	44 g



Alle Angaben ohne Gewähr.
All information is provided without guarantee.

IPL safety eyewear ONTOR Shade 5 (suitable also for spectacle wearers) offer safety for operator and patient during IPL treatments.

Properties:

Our **IPL safety goggles ONTOR** are lightweight, versatile **safety goggles** made of a plastic frame, which is perfectly suitable as over-glasses for spectacle wearers. The frame offers a very large field of vision and wide side coverage, with excellent wearing comfort due to flexible temples adjustable in length and inclination, and a special soft-rubber padding in the forehead area. **IPL protective goggles ONTOR** provide safety for users and patients during **IPL treatments**.

- Attenuation: **Shade 5** according **EN 169**
- suitable for **IPL applications**
- modern sport design
- high wearing acceptance due to perfect fit
- unrestricted field of vision
- suitable for spectacle wearers
- Individually adjustable temple length
- Adjustable inclination
- lightweight (44g)

Included in delivery: microfibre bags with PROTECT-logo, polybag, user's manual, cleaning cloth with PROTECT-logo, glasses head strap