


IPL protection eyewear ONTOR F-4

Product information:

Item no.:	600-K234-ONTO-54
Filter color:	grey
Frame color:	black/white
Filter material:	special plastic
Product weight:	44 g



Alle Angaben ohne Gewähr.
All information is provided without guarantee.

IPL safety eyewear ONTOR F4 (suitable also for spectacle wearers) offer safety for operator and patient during IPL treatments.

Properties:

Our **IPL safety goggles ONTOR** are lightweight, versatile **safety goggles** made of a plastic frame, which is perfectly suitable as over-glasses for spectacle wearers. The frame offers a very large field of vision and wide side coverage, with excellent wearing comfort due to flexible temples adjustable in length and inclination, and a special soft-rubber padding in the forehead area. **IPL protective goggles ONTOR** provide safety for users and patients during **IPL treatments**.

- suitable for **IPL applications**
- very large field of vision

Certified according to **EN 12609-1**

- Protection level: **F-4**

Maximum spectral transmittance in **UV** and near UV range:

- 250 to 315 nm: 0.1
- >315 to 380 nm: 0.4 %
- >380 to 450 nm: $\tau(\text{VIS}, \text{D65})$

Light transmittance $\tau(\text{VIS}, \text{D65})$:

- $3,2\% \leq \tau(\text{VIS}, \text{D65}) < 8,5\%$

Maximum spectral transmittance in the **IR** range:

- 780 nm to 3000 nm: 50

Included in delivery: microfibre bags with PROTECT-logo, polybag, user's manual, cleaning cloth with PROTECT-logo, glasses head strap