

## Eyecaps metal with offset handle M

### Product information:

**Item no.:** 600-AUGK-M-APP-M  
**Product weight:** 0 g

Application for all laser wavelengths Material:  
Stainless steel, opaque Dimensions: approx. 25.5 x  
23.0 mm One pair with fixed handle Eye caps  
autoclavable



### Properties:

Application for all laser wavelengths  
Material: Stainless steel, opaque  
Dimensions: approx. 25.5 x 23.0 mm  
One pair with fixed handle  
Eye caps autoclavable

Similar to contact lenses, the eye caps are placed directly on the eye. To prevent injury, the eye caps have rounded edges with a highly polished interior and a matte finish on the exterior to protect medical personnel from reflections.

The stainless steel variant is 100% opaque, making it suitable protection for all laser wavelengths.

Other variants, such as [plastic models](#), are also available.

**EYE CAPS ARE NOT DISCOUNTABLE!**